



# Guía de selección de desinfectante para el coronavirus

Nombre & Número	Ingrediente Activo	Dilución	pH	Acinetobacter	Enterobacteriaceae (carbapenem-resistance)	Hepatitis A	Hepatitis B	Hepatitis C	HIV	Human Coronavirus	Human Coronavirus (SARS Associated)	Influenza	Klebsiella	MRSA	Norovirus	Pseudomonas aeruginosa	Staphylococcus aureus	VRE	Tuberculosis (TB)	VISA	VRSA	Canine Parvo	Salmonella enterica/choleeraesuis
Triforce (333)	Quat	1:256	Alkaline	3 min			3 min	3 min	1 min	3 min	10 min	1 min	3 min	3 min		3 min	3 min	3 min		3 min			3 min
Quat-Stat 5™ (341)	Quat	1:256	Alkaline	5 min	5 min		5 min	5 min	5 min	5 min	5 min	5 min	5 min	5 min	5 min	5 min	5 min	5 min			5 min		5 min
Fight-Bac® RTU (311)	Quat	RTU	Alkaline			10 min	1 min	1 min	1 min	2 min	2 min	2 min		3 min	30 sec	3 min	3 min	3 min	5 min	3 min		10 min	3 min
pH7Q (316)	Quat	1:64	Neutral	10 min			10 min	10 min	4 min	1 min	10 min	2 min	10 min	10 min		10 min	10 min	10 min		10 min			10 min
pH7Q Dual (355)	Quat	1:256	Neutral	10 min	10 min		10 min	10 min	10 min	10 min	10 min	10 min	10 min	10 min		10 min	10 min	10 min		10 min	10 min	10 min	10 min
Pine Quat (304)	Quat	1:32	Neutral	10 min			10 min	10 min	4 min	10 min	10 min	10 min	10 min	10 min			10 min	10 min		10 min			10 min
Symlicity™ Sanibet™ Multi-Range (237)	Quat	3 oz/5 gal	Neutral				10 min	10 min	10 min	10 min	10 min	10 min	10 min	10 min	10 min	10 min	10 min						10 min

- Triforce Sanitizer for soft, non-food contact surfaces at 1:256